

Product Datasheet

Human ENA-78/CXCL5 ELISA Kit (orb339793)

Description	Human Epithelial neutrophil activating peptide 78 (ENA-78/CXCL5) ELISA Kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	CXCL5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P42830
Sensitivity	0.567 ng/mL
Expiration Date	Please enquire.