

Product Datasheet

Human GM-CSF ELISA Kit (orb339796)

Description	Human Granulocyte-Macrophage Colony Stimulating Factor, GM-CSF ELISA Kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	CSF2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P04141
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.