

Product Datasheet

Human IL-8 ELISA Kit (orb339804)

Description	Human Interleukin 8, IL-8 ELISA KIT
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	IL8
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, cell culture supernates, saliva, urine, cerebrospinal fluid (CSF), tissue homogenates, cell lysates
Uniprot ID	P10145
Sensitivity	7.110 pg/mL
Expiration Date	Please enquire.