

Product Datasheet

Human MDC/CCL22 ELISA Kit (orb339809)

Description	Human Macrophage-Derived Chemokine (MDC/CCL22) ELISA kit
Reactivity	Human
Range	15.625 pg/mL-1000 pg/mL
Target	CCL22
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	000626
Sensitivity	10.324 pg/mL
Expiration Date	Please enquire.