

Product Datasheet

Human sCD14 ELISA Kit (orb339816)

Description	Human soluble cluster of differentiation 14, sCD14 ELISA Kit
Reactivity	Human
Range	0.234 ng/mL-15 ng/mL
Target	CD14
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, urine
Uniprot ID	P08571
Sensitivity	0.229 ng/mL
Expiration Date	Please enquire.