

Product Datasheet

Human TNFRSF14 ELISA Kit (orb339818)

Description	Human Tumor necrosis factor receptor superfamily member 14 (TNFRSF14) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	TNFRSF14
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q92956
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.