

## Product Datasheet

### Human TNFRSF17 ELISA Kit (orb339819)

<b>Description</b>	Human Tumor necrosis factor receptor superfamily member 17 (TNFRSF17) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	7.8 pg/ml- 500 pg/ml
<b>Target</b>	TNFRSF17
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q02223</a>
<b>Sensitivity</b>	1.95 pg/ml
<b>Expiration Date</b>	Please enquire.