

Product Datasheet

Mouse Flt-3L ELISA Kit (orb339824)

Description	Mouse FMS-like tyrosine kinase 3 ligand, Flt-3L ELISA Kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	FLT3LG
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P49772
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.