

Product Datasheet

Mouse MDC ELISA Kit (orb339827)

Description	Mouse Macrophage-Derived Chemokine, MDC ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	CCL22
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates, cell lysates
Uniprot ID	O88430
Sensitivity	23.672 pg/mL
Expiration Date	Please enquire.