

Product Datasheet

Mouse MCP-5 ELISA Kit (orb339828)

Description Mouse Monocyte Chemotactic Protein 5 (MCP-5) ELISA Kit

Reactivity Mouse

Range 15.6 pg/mL-1000 pg/mL

Target MCP-5

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma, tissue homogenates

Uniprot ID [Q62401](#)

Sensitivity 16.82 pg/mL

Expiration Date Please enquire.