

Product Datasheet

Pig TNF/TNFA/TNFSF2 ELISA Kit (orb339833)

Description	Pig Tumor necrosis factor (TNF/TNFA/TNFSF2) ELISA kit
Reactivity	Porcine
Range	0.047 ng/mL-3 ng/mL
Target	TNF
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P23563
Sensitivity	0.040 ng/mL
Expiration Date	Please enquire.