

Product Datasheet

Bovine LGB ELISA Kit (orb339837)

Description	Bovine beta-lactoglobulin (Beta-LG) (allergen Bos d 5) (LGB) ELISA kit
Reactivity	Bovine
Range	0.156 µg/mL-10 µg/mL
Target	LGB
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, milk
Uniprot ID	P02754
Sensitivity	0.129 µg/mL
Expiration Date	Please enquire.