

## Product Datasheet

### Guinea Pig Alb ELISA Kit (orb339853)

<b>Description</b>	Guinea Pig Albumin (Alb) ELISA kit
<b>Reactivity</b>	Guinea pig
<b>Range</b>	0.156 µg/mL-10 µg/mL
<b>Target</b>	ALB
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q6WDN9</a>
<b>Sensitivity</b>	0.09 µg/mL
<b>Expiration Date</b>	Please enquire.