

## Product Datasheet

### Human NGAL ELISA Kit (orb339868)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Human neutrophil gelatinase-associated lipocalin, NGAL ELISA Kit |
| <b>Reactivity</b>        | Human  |
| <b>Range</b>             | 15.6 pg/mL-1000 pg/mL  |
| <b>Target</b>            | NGAL   |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                     |
| <b>Sample Types</b>      | serum, plasma, urine, tissue homogenates                         |
| <b>Uniprot ID</b>        | <b>P80188</b>  |
| <b>Sensitivity</b>       | 3.9 pg/mL  |
| <b>Expiration Date</b>   | Please enquire.  |