

Product Datasheet

Rat MAU/ALB ELISA Kit (orb339896)

| | |
|--------------------------|---|
| Description | Rat microalbuminuria (MAU/ALB) ELISA kit |
| Reactivity | Rat |
| Range | 0.01 µg/mL-40 µg/mL |
| Target | MAU/ALB |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, urine |
| Uniprot ID | P02770 |
| Sensitivity | 0.01 µg/mL |
| Expiration Date | Please enquire. |