

## Product Datasheet

### Mouse HBeAb ELISA Kit (orb339986)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Mouse hepatitis B virus e antibody, HBeAb ELISA Kit |
| <b>Reactivity</b>        | Mouse   |
| <b>Range</b>             | Request Information                                 |
| <b>Target</b>            | HBeAb   |
| <b>Note</b>              | For research use only                               |
| <b>Application notes</b> | Detection Wavelength: 450 nm                        |
| <b>Sample Types</b>      | serum, plasma                                       |
| <b>Sensitivity</b>       | Request Information                                 |
| <b>Expiration Date</b>   | Please enquire.                                     |