

## Product Datasheet

### Human CSTB/CST6/STFB ELISA Kit (orb340024)

<b>Description</b>	Human Cystatin-B (CSTB/CST6/STFB) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	0.313 ng/mL-20 ng/mL
<b>Target</b>	CSTB
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, urine, saliva, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<b>P04080</b>
<b>Sensitivity</b>	0.233 ng/mL
<b>Expiration Date</b>	Please enquire.