

## Product Datasheet

### Human NOV ELISA Kit (orb340030)

<b>Description</b>	Human Protein NOV homolog (NOV) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	62.5 pg/mL-4000 pg/mL
<b>Target</b>	NOV
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P48745</b>
<b>Sensitivity</b>	15.6 pg/mL
<b>Expiration Date</b>	Please enquire.