

## Product Datasheet

### Canine Testosterone ELISA Kit (orb340052)

<b>Description</b>	Canine Testosterone (T) ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	0.1 ng/mL-20 ng/mL
<b>Target</b>	T
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Sensitivity</b>	0.05 ng/mL
<b>Expiration Date</b>	Please enquire.