

## Product Datasheet

# Human Chitinase-3-like Protein 1 ELISA Kit (orb340082)

<b>Description</b>	Human Chitinase-3-like Protein 1 (YKL-40/CHI3L1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	46.875 pg/mL-3000 pg/mL
<b>Target</b>	CHI3L1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">P36222</a>
<b>Sensitivity</b>	78.065 pg/mL
<b>Expiration Date</b>	Please enquire.