

Product Datasheet

Human Decorin/Bone proteoglycan II ELISA Kit (orb340083)

Description	Human Decorin/Bone proteoglycan II (DCN) ELISA Kit
Reactivity	Human
Range	3.9 pg/mL-250 pg/mL
Target	DCN
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates.
Uniprot ID	P07585
Sensitivity	0.97 pg/mL
Expiration Date	Please enquire.