

Product Datasheet

Human VEGFR-3/Flt-4 ELISA Kit (orb340100)

Description	Human Vasclar endothelial cell growth factor receptor 3, VEGFR-3/Flt-4 ELISA Kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	VEGFR-3/Flt-4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P35916
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.