

Product Datasheet

Mouse Decorin ELISA Kit (orb340107)

| | |
|--------------------------|-----------------------------------|
| Description | Mouse Decorin (DCN) ELISA kit |
| Reactivity | Mouse |
| Range | 78 pg/mL-5000 pg/mL |
| Target | DCN |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | P28654 |
| Sensitivity | 19.5 pg/mL |
| Expiration Date | Please enquire. |