

## Product Datasheet

### Rat Estradiol ELISA Kit (orb340124)

<b>Description</b>	Rat Estradiol, E2 ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	50 pg/ml-1200 pg/ml.
<b>Target</b>	E2
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma
<b>Sensitivity</b>	40 pg/ml.
<b>Expiration Date</b>	Please enquire.