

Product Datasheet

Rat IGF-1 ELISA Kit (orb340125)

Description	Rat Insulin-like growth factor 1, IGF-1 ELISA Kit
Reactivity	Rat
Range	0.156 ng/mL-10 ng/mL
Target	IGF1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P08025
Sensitivity	0.151 ng/mL
Expiration Date	Please enquire.