

## Product Datasheet

### Rat PROG ELISA Kit (orb340126)

<b>Description</b>	Rat Progesterone, PROG ELISA kit
<b>Reactivity</b>	Rat
<b>Range</b>	1 ng/ml-70 ng/ml
<b>Target</b>	PROG
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates.
<b>Sensitivity</b>	0.2 ng/ml
<b>Expiration Date</b>	Please enquire.