

Product Datasheet

Rat Testosterone ELISA Kit (orb340127)

| | |
|--------------------------|---|
| Description | Rat Testosterone, T ELISA Kit |
| Reactivity | Rat |
| Range | 0.13 ng/mL-25.6 ng/mL |
| Target | T |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, cell culture supernates, tissue homogenates |
| Sensitivity | 0.06 ng/mL |
| Expiration Date | Please enquire. |