

## Product Datasheet

### Sheep FSH ELISA Kit (orb340131)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Sheep Follicle Stimulating Hormone (FSH) ELISA kit |
| <b>Reactivity</b>        | Sheep  |
| <b>Range</b>             | Request Information                                |
| <b>Target</b>            | FSH  |
| <b>Note</b>              | For research use only                              |
| <b>Application notes</b> | Detection Wavelength: 450 nm                       |
| <b>Uniprot ID</b>        | <a href="#">P01227</a>                             |
| <b>Sensitivity</b>       | Request Information                                |
| <b>Expiration Date</b>   | Please enquire.                                    |