

## Product Datasheet

### Sheep Luteinizing Hormone ELISA Kit (orb340132)

<b>Description</b>	Sheep Luteinizing Hormone (LH) ELISA Kit
<b>Reactivity</b>	Sheep
<b>Range</b>	Request Information
<b>Target</b>	LH
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P01231</b>
<b>Sensitivity</b>	Request Information
<b>Expiration Date</b>	Please enquire.