

Product Datasheet

Lestaurtinib (CEP701) (orb340482)

Description	Lestaurtinib (CEP701)
Target	Kinase
Form/Appearance	powder
Storage	-20°C
Note	For research use only
Formula	$C_{26}H_{21}N_3O_4$
Purity	> 98%
MW	439.46
CAS Number	111358-88-4
Solubility (25°C)	Soluble in DMSO (up to 100 mg/ml) or in Ethanol (up to 20 mg/ml)
Expiration Date	12 months from date of receipt.