

## Product Datasheet

### Sorafenib (orb340489)

<b>Description</b>	Sorafenib
<b>Target</b>	Kinase
<b>Form/Appearance</b>	powder
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Formula</b>	$C_{21}H_{16}ClF_3N_4O_3$
<b>Purity</b>	> 98%
<b>MW</b>	464.83
<b>CAS Number</b>	284461-73-0
<b>Solubility (25°C)</b>	Soluble in DMSO (up to 200 mg/ml) or in Ethanol (up to 3 mg/ml).
<b>Expiration Date</b>	12 months from date of receipt.