



Product Datasheet

FBL antibody (orb340733)

Description	Rabbit polyclonal antibody to FBL
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Recombinant full length protein of human Fibrillarin
Target	FBL
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P22087 , P22509 , P35550
Entrez	2091 , 292747 , 102643269 , 14113
Dilution Range	WB: 1:500-1:2000, IF/ICC: 1:50-1:200, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of Fibrillarin exp...



Immunofluorescent analysis of Fibrillari...