



Product Datasheet

TUBA4A antibody (orb340868)


Description	Rabbit polyclonal antibody to TUBA4A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant full length protein of human Alpha-tubulin 4a
Target	TUBA4A
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P68366, P68368, Q5XIF6
Entrez	22145, 7277, 316531
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:25-1:50, IF/ICC: 1:25-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of Alpha-tubulin 4a



Immunohistochemical analysis of Alpha-tubulin 4a



Immunofluorescent analysis of Alpha-tubulin 4a