

## Product Datasheet

### PKC-412 (orb341511)

<b>Description</b>	PKC-412 is a broad spectrum protein kinase inhibitor. Inhibits conventional PKC isoforms ( $\alpha$ , $\beta$ , $\gamma$ ...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{35}H_{30}N_4O_4$
<b>MW</b>	570.6371
<b>CAS Number</b>	[120685-11-2]
<b>SMILES</b>	<chem>C[C@@]12[C@@H]([C@@H](C[C@@H](O1)N3C4=CC=CC=C4C5=C6C(=C7C8=CC=CC=C8N2C7=C53)CNC6=O)N(C)C(=O)C9=CC=CC=C9)OC</chem>
<b>Expiration Date</b>	12 months from date of receipt.