


Product Datasheet

AB-MECA (orb341594)

Description	AB-MECA is a high affinity A3 adenosine receptor agonist, has high affinity for recombinant A1 and A3 receptors.
Note	For research use only
Formula	$C_{18}H_{21}N_7O_4$
MW	399.403842687607
CAS Number	[152918-26-8]
SMILES	<chem>NC1C=CC(CNC2=NC=NC3N([C@@H]4O[C@H](C(NC)=O)[C@@H](O)[C@H]4O)C=NC2=3)=CC=1</chem>
Expiration Date	12 months from date of receipt.



Chemical structure of AB-MECA

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)