


## Product Datasheet

### Sal003 (orb341656)

<b>Description</b>	Sal003 is a potent cell-permeable analog of the eIF2 $\alpha$ phosphatase inhibitor Salubrinal with enhance...
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>18</sub> H <sub>15</sub> Cl <sub>4</sub> N <sub>3</sub> OS
<b>MW</b>	463.21
<b>CAS Number</b>	[1164470-53-4]
<b>SMILES</b>	O=C(C=CC1C=CC=CC=1)NC(C(Cl)(Cl)Cl)NC(=S)NC1C=CC(Cl)=CC=1
<b>Expiration Date</b>	12 months from date of receipt.



Chemical structure of Sal003