


## Product Datasheet

### NVP-BEZ235 Tosylate (orb341702)

<b>Description</b>	BEZ235 (NVP-BEZ235) is a dual ATP-competitive PI3K and mTOR inhibitor of p110 $\alpha$ , p110 $\gamma$ , p110 $\delta$ and ...
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>30</sub> H <sub>23</sub> N <sub>5</sub> O <sub>3</sub> S
<b>MW</b>	641.73814
<b>CAS Number</b>	[1028385-32-1]
<b>SMILES</b>	<chem>CN(C1=C2C3=CC(C4=CC5=CC=CC=C5N=C4)=CC=C3N=C1)C(N2C6=CC=C(C=C6)C(C)(C#N)C)=O.CC7=CC=C(S(=O)(O)=O)C=C7</chem>
<b>Expiration Date</b>	12 months from date of receipt.



Chemical structure of NVP-BEZ235 Tosylat...