

## Product Datasheet

### Pimavanserin(ACP-103) (orb341711)

<b>Description</b>	Pimavanserin(ACP-103) is a potent and selective 5-HT <sub>2A</sub> receptor inverse agonist with mean pIC <sub>50</sub> of w...
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>25</sub> H <sub>34</sub> FN <sub>3</sub> O <sub>2</sub>
<b>MW</b>	427.56
<b>CAS Number</b>	[706779-91-1]
<b>SMILES</b>	<chem>N(CC1=CC=C(F)C=C1)(C1CCN(C)CC1)C(NCC1=CC=C(OCC(C)C)C=C1)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.