

## Product Datasheet

### SKF38393 HCl (orb341772)

<b>Description</b>	SKF38393 HCl is a selective dopamine D1/D5 receptor agonist.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{16}H_{17}NO_2 \cdot HCl$
<b>MW</b>	291.77
<b>CAS Number</b>	[62717-42-4]
<b>SMILES</b>	<chem>Cl.C1(C2C3C=C(O)C(O)=CC=3CCNC2)C=CC=CC=1</chem>
<b>Expiration Date</b>	12 months from date of receipt.