



Product Datasheet

Epitope Tag Antibody Sampler Kit (orb342265)


Description	Epitope Tag Antibody Sampler Kit
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, FC, IF, IHC, IP, WB
Immunogen	GFP, GST, DYKDDDDK peptide (the same epitope recognized by Sigma's Anti-FLAG antibodies), 6X His peptide
Storage	Store Epitope Tag Kit at -20° C prior to opening. Aliquot antibodies and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	Epitope Tag kit contains: 6X HIS Antibody 100 µg DYKDDDDK Antibody (Binds to same epitope as Sigma's Anti-FLAG M2 Antibody) 250 µg GFP Antibody 100 µg GST Antibody 1.0mg Anti-RABBIT IgG (H&L) (GOAT) Antibody Peroxidase Conjugated 100 µg Anti-GOAT IgG (H&L) (RABBIT) Antibody Peroxidase Conjugated 100 µg
Purity	6X HIS antibody, GFP Antibody, DYKDDDDK Antibody (Binds to same epitope as Sigma's Anti-FLAG M2 Antibody), GST Antibody and Conjugated Secondary Antibodies are all affinity purified.
Dilution Range	ELISA: User Optimized, ChIP: User Optimized, FC: User Optimized, IHC: User Optimized, IF: User Optimized, IP: User Optimized, WB: User Optimized
Expiration Date	Please enquire.



Anti-6X His epitope tag polyclonal antibody...



Biorbyt Affinity Purified anti FLAG™ P...



Biorbyt anti GFP polyclonal antibody (or...