

Product Datasheet

Western Blot Kit (orb342271)

Description Western Blot Kit

Reactivity Mouse

Conjugation Unconjugated

Tested Applications WB

Preservatives Wash buffers MUST NOT contain SODIUM AZIDE or other inhibitors of peroxidase activity!

Form/Appearance n/a

Concentration n/a

Storage See kit insert for complete instructions.


Note For research use only

Application notes Biorbyt' Chemiluminescent Western Blot Kit, namely our FemtoMax Chemiluminescent Substrate Kit for Western Blotting, combines all of the necessary reagents with a rapid proven protocol and the extremely high signal detection of our luminol substrate. The Chemiluminescent Western Blot Kit design includes straightforward procedures and color-coding to allow for ease of use. This kit contains sufficient substrate for up to 30 mini blots at 7.5 x 8 cm² and is stable for at least 1 year when stored at +4°C.

Purity This kit contains: anti-mouse peroxidase, mouse IgG control, BSA, tween-20, and Femtomax, along with the protocol and additional supplies. See kit insert for complete details and instructions.

Dilution Range WB: 1:20,000

Expiration Date Please enquire.



Chemiluminescent Western Blot Kit shows ...

