


## Product Datasheet

### Western Blot Kit (orb342272)

<b>Description</b>	Western Blot Kit
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Preservatives</b>	Wash buffers MUST NOT contain SODIUM AZIDE or other inhibitors of peroxidase activity!
<b>Form/Appearance</b>	n/a
<b>Storage</b>	See kit insert for complete instructions.
<b>Note</b>	For research use only
<b>Application notes</b>	Biorbyt' Chemiluminescent Western Blot Kit, namely our FemtoMax Chemiluminescent Substrate Kit for Western Blotting, combines all of the necessary reagents with a rapid proven protocol and the extremely high signal detection of our luminol substrate. The Chemiluminescent Western Blot Kit design includes straightforward procedures and color-coding to allow for ease of use. This kit contains sufficient substrate for up to 30 mini blots at 7.5 x 8 cm <sup>2</sup> and is stable for at least 1 year when stored at +4°C.
<b>Purity</b>	This kit contains: anti-rabbit peroxidase, rabbit IgG control, BSA, tween-20, and Femtomax, along with the protocol and additional supplies. See kit insert for complete details and instructions.
<b>Dilution Range</b>	WB: User Optimized
<b>Expiration Date</b>	Please enquire.



FemtoMax™  
Chemiluminescent  
Western Blo...

