

Product Datasheet

Biotin anti-Goat IgG Antibody MaxTag Histo for IHC (orb342291)

Description	Biotin anti-Goat IgG Antibody MaxTag Histo for IHC
Conjugation	Biotin
Tested Applications	IHC
Immunogen	Goat IgG whole molecule
Preservatives	0.01% (w/v) Sodium Azide and 0.01% (w/v) Gentamicin Sulfate
Form/Appearance	Liquid (sterile filtered)
Storage	See kit insert for complete instructions.
Note	For research use only
Expiration Date	12 months from date of receipt.