

## Product Datasheet

# Biotin Western Blot Anti-Goat IgG Antibody Maxtag (orb342304)

<b>Description</b>	Biotin Western Blot Anti-Goat IgG Antibody Maxtag
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Goat IgG whole molecule
<b>Preservatives</b>	0.01% (w/v) Sodium Azide and 0.01% (w/v) Gentamicin Sulfate
<b>Form/Appearance</b>	Liquid (sterile filtered)
<b>Concentration</b>	1.0 mg/mL
<b>Storage</b>	See kit insert for complete instructions.
<b>Note</b>	For research use only
<b>Application notes</b>	See Kit insert for complete instructions.
<b>Expiration Date</b>	12 months from date of receipt.