



Product Datasheet

Cyclin Antibody Sampler Kit (orb342331)


Description	Cyclin Antibody Sampler Kit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Form/Appearance	Liquid (sterile filtered)
Concentration	n/a
Storage	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	Each kit contains one unit of the following products: (100 µL) Anti-Cyclin A (RABBIT) Antibody (100 µL) Anti-Cyclin B1 (RABBIT) Antibody (100 µL) Anti-Cyclin D1 (RABBIT) Antibody (100 µL) Anti-Cyclin E (RABBIT) Antibody (100 µg) Anti-RABBIT IgG (H&L) (GOAT) Antibody Peroxidase Conjugated
Dilution Range	ELISA: 1:2,000 - 1:10,000, IP: 1:100, WB: 1:500 - 1:5,000
Expiration Date	Please enquire.



Biorbyt's Anti-Cyclin A antibody was dil...



Biorbyt's anti-Cyclin B1 antibody was di...



Western blot analysis is shown using Bio...