

## Product Datasheet

### GLB1L2 antibody (HRP) (orb34361)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | GLB1L2 antibody (HRP)   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human   |
| <b>Conjugation</b>         | HRP   |
| <b>Tested Applications</b> | ELISA   |
| <b>Immunogen</b>           | Recombinant Human Beta-galactosidase-1-like protein 2 protein (401-636AA) |
| <b>Target</b>              | GLB1L2  |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol            |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.             |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Uniprot ID</b>          | <a href="#">Q8IW92</a>  |
| <b>Dilution Range</b>      | ELISA   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |