

Product Datasheet

NFKB p65 (Rel A) (phospho-S276) peptide (orb343845)

| | |
|--------------------------|---|
| Description | NFKB p65 (Rel A) (phospho-S276) peptide |
| Conjugation | Unconjugated |
| Preservatives | 0.01% (w/v) Sodium Azide |
| Form/Appearance | Liquid (sterile filtered) |
| Concentration | 1.0 mg/mL |
| Storage | Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. |
| Note | For research use only |
| Application notes | Control peptide should be used at 1.0 µg per 1.0 µl of antiserum in per assay. |
| Purity | Greater than 95% specific peptide |
| Source | Human |
| Uniprot ID | Q04206 |
| NCBI | 223468676 |
| Expiration Date | 12 months from date of receipt. |