



Product Datasheet

H3 (phospho-Ser10) antibody (orb34388)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	H3C1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	15404
Uniprot ID	P68431
NCBI	NP_003520.1 , NP_003524.1 , NP_003521.2 , NP_003526.1 , NP_066298.1 , NP_003523.1 , NP_003522.1 , NP_003528.1 , NP_003525.1 , NP_003527.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of in CEM cell lin...