

## Product Datasheet

### Apg7 peptide (orb343880)

<b>Description</b>	Apg7 peptide
<b>Conjugation</b>	Unconjugated
<b>Preservatives</b>	None
<b>Form/Appearance</b>	Liquid (sterile filtered)
<b>Concentration</b>	1.0 mg/mL
<b>Storage</b>	Store vial at -20 °C prior to opening. Aliquot contents and freeze at -20 °C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4 °C as an undiluted liquid. Dilute only prior to immediate use.
<b>Note</b>	For research use only
<b>Application notes</b>	Apg7 Control peptide should be used at 1.0 µg per 1.0 µl of antiserum in per assay.
<b>Purity</b>	Greater than 95% specific peptide
<b>Source</b>	Human
<b>Uniprot ID</b>	<a href="#">O95352</a>
<b>Expiration Date</b>	12 months from date of receipt.