



## Product Datasheet

### H3 (phospho-Ser10) antibody (orb34389)

|                            |                                                                                                                                                                                                                                                                                                           |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Peptide affinity purified rabbit polyclonal antibody.                                                                                                                                                                                                                                                     |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                                                                                                    |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                                                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                                                                              |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                                                                                                                                                                                                                 |
| <b>Target</b>              | H3C1                                                                                                                                                                                                                                                                                                      |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.                                                                                                                              |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles                                                                                                                                                                    |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                     |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                                                                                                                                                |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                                                                |
| <b>MW</b>                  | 15404                                                                                                                                                                                                                                                                                                     |
| <b>Uniprot ID</b>          | <a href="#">P68431</a>                                                                                                                                                                                                                                                                                    |
| <b>NCBI</b>                | <a href="#">NP_003520.1</a> , <a href="#">NP_003524.1</a> , <a href="#">NP_003521.2</a> , <a href="#">NP_003526.1</a> , <a href="#">NP_066298.1</a> , <a href="#">NP_003523.1</a> , <a href="#">NP_003522.1</a> , <a href="#">NP_003528.1</a> , <a href="#">NP_003525.1</a> , <a href="#">NP_003527.1</a> |
| <b>Dilution Range</b>      | WB: 1:1000, IHC-P: 1:50-100                                                                                                                                                                                                                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                                                                           |



Western blot analysis of in HepG2 cell l...